

BIRR and IVF Labware: it's embryologic!

smART 5-Well Dish

5 good reasons to choose BIRR's 5-Well



- 1 Easier handling and access*
- 2 Improved visualization**
- 3 Large labeling area***
- 4 Consistent heat transfer*
- 5 MEA+ tested


* The well's larger diameter reduces the angle at which the pipette needs to be held during specimen manipulation, facilitating operations.

** The bottom of the dish is raised by 0.2mm, preventing scratches on the bottom of the dish and minimizing potential negative effects of temperature surges from the heated surface. Round edges at the bottom of the well minimise shading and improve specimen visualisation.

*** Large labelling area consents clear labelling with patients details and RFID tags without obstructing specimen visualization.

Want to learn more about the smART characteristics of our 5-Well Dish?

Please scan this QR code!



Our smART Labware introduces innovative features and added value on ease of use

BIRR It's embryologic!