

BIRR and IVF Labware: it's embryologic!

Why choose the 4+8-Well Dish



Two different drop sizes to accommodate flexible usage for: cell culture (gametes or embryos), oocyte insemination (IVF), cumulus-oocyte complex denudation, handling or manipulation of gametes.




- 1 New 'ears' design**
 - Provides larger labeling space
 - Enhances grip for improved handling
- 2 Concave well design**
 - Simplifies embryo visibility
- 3 Raised well bottom**
 - Prevents scratching
 - Buffers temperature changes
- 4 Well-numbering 1-12**
 - Supports tracking of embryos

Our smART Labware introduces innovative features and added value on ease of use

smART
LABWARE

SCAN to learn
more about the
BIRR dishes portfolio



BIRR It's embryologic!