

Serological Pipettes


Guidelines For Use

Serological pipettes offer effectively improve experimental efficiency when used in fume hoods, biosafety cabinets, isolators and other confined spaces.



PakGent® 's rigorous quality specifications and convenient functionality make PakGent® serological pipettes an excellent choice for routine laboratory work and cell studies. Bold, clear bidirectional scale makes liquid capacity and sample addition easy to read. A variety of specifications and packaging options are available. Certified scale accuracy and precision calibration ensures +/-2% accuracy in the 1mL and 2mL total volume range and +/-1% accuracy in the 5 mL, 10mL, 25mL, 50mL total volume range.

Features	Description
Material	USP Class VI polystyrene
Sterility	SAL of 10 ⁻⁶
Pack	Optional Individually packed pipettes



● 1ml Serological Pipette, PS, Sterile

Part	Package Type	Pack
SPTS01	Individual pack	50/bag, 20bags/case
SPTB01	bulk pack	50/bag, 20bags/case

● 5ml Serological Pipette, PS, Sterile

Part	Package Type	Pack
SPTS05	Individual pack	50/bag, 10bags/case
SPTB05	bulk pack	50/bag, 10bags/case

● 10ml Serological Pipette, PS, Sterile

Part No.	Pcs	Package Type	Pack
SPTS10	one-pc	Individual pack	50/bag, 10bags/case
SPTB10	one-pc	bulk pack	50/bag, 10bags/case
SPTS10A	two-pcs	Individual pack	50/bag, 10bags/case
SPTB10A	two-pcs	bulk pack	50/bag, 10bags/case

● 2ml Serological Pipette, PS, Sterile



Part	Package Type	Pack
SPTS02	Individual pack	50/bag, 20bags/case
SPTB02	bulk pack	50/bag, 20bags/case

● 25ml Serological Pipette, PS, Sterile

Part	Package Type	Pack
SPTS25	Individual pack	25/bag, 10bags/case
SPTB25	bulk pack	25/bag, 10bags/case

● 50ml Serological Pipette, PS, Sterile

Part	Package Type	Pack
SPTS50	Individual pack	25/bag, 8bags/case
SPTB50	bulk pack	25/bag, 8bags/case



SIZE (mL)	Appearance size				Pipe tip			Pipe tail			Extrusion pipe		
	A	B	C	D	L1	ID1	OD1	L3	ID3	OD3	L2	ID2	OD2
1	266±2	2.3±0.3	4.6±0.2	4.6±0.2	-	1.6±0.3	2.3±0.3	-	2.8±0.2	4.6±0.2	-	2.8±0.2	4.6±0.2
2	270±2	2.7±0.3	6.3±0.2	6.3±0.2	-	1.7±0.3	2.7±0.3	-	4.1±0.2	6.3±0.2	-	4.1±0.2	6.3±0.2
5	342±2	3.2±0.3	8.1±0.2	8.1±0.2	-	2.2±0.3	3.2±0.3	-	6±0.2	8.1±0.2	-	6±0.2	8.1±0.2
10	344±2	3.6±0.15	7.9±0.2	10±0.2	37	1.5±0.1	3.6±0.15	19.2±0.2	5.5±0.2	7.9±0.2	285	8.2±0.2	10±0.2
25	313±2	5±0.2	8.1±0.2	15±0.2	35	2.5±0.2	5±0.2	21.7±0.2	5.2±0.2	8.1±0.2	255	13±0.2	15±0.2
50	353±2	5.6±0.2	8.1±0.2	18±0.2	35	3±0.2	5.6±0.2	21.7±0.2	5.2±0.2	8.1±0.2	295	16±0.2	18±0.2

PakGent®

PakGent Bioscience (Suzhou) Co., Ltd.

Address: No. 9 Lecheng Road, Bixi Street, Changshu, Jiangsu Province, China

Tel: +86 (0) 512 6670 0017

Email: info@pakgentbio.com

Website: www.pakgentbio.com



Catalog