

1000 µl Blue Graduated TipOne Tip, Racked (ns)

TipOne® is the original refillable, reusable and recyclable pipette tip system. Recent improvements to the design of the system, based on customer-feedback, is why TipOne® is one of the leading tips for research labs.

Convenience and quality are at the heart of the system and TipOne® remains the same tip that you have used and trusted over the years.

Tip fit

One tip fits all! The tip collar is engineered to fit all popular pipettes.

Manufacture

The use of multiple cavity tools for the injection moulding process guarantees high conformity within one batch. The automated moulding and assembling procedures are carried out under clean-room conditions to eliminate contamination.

Unique collar

Our newest tips can be recognised by the jagged design on the tip collar.

Material

All TipOne® tips are manufactured from pure, virgin polypropylene to ensure the highest quality product.

Green tip design

A reduction of up to 40 % in the material used for the production of TipOne® tips helps to save valuable resources and is beneficial for the environment.

Highly polished moulds

For an extremely smooth surface to reduce sample retention.


Quality Control

Automated Quality Control procedures ensure that TipOne® is inspected for 100% performance. Every TipOne® product is certified RNase, DNase, DNA and Pyrogen free.

Product Features

- **Graduated** at 250 µl, 500 µl and 1000 µl
- **Universal fit:** One tip fits all! The tip collar is engineered to fit all popular pipettes.
- **Highly polished moulds:** For an extremely smooth surface to reduce sample retention.
- **Material:** All TipOne® tips are manufactured from pure, virgin polypropylene to ensure the highest quality product.
- **Certified:** Every TipOne® product is certified RNase, DNase, DNA and Pyrogen free.



1000 µl Blue Graduated TipOne Tip, Racked (ns)			
	Description	Pack Size	Cat. No.
	1000 µl Blue Graduated TipOne Tip, Racked (ns)	960 Tips (10 Racks × 96 Tips)	S1111-6801

Technical Specifications	
Autoclavable	Yes
Color Code Wafer	Blue
Material	Polypropylene
Pack Size	960 Tips (10 Racks × 96 Tips)
Sterile	No
Color	Blue
Volume	1,000 µl
For use with	Abimed®, Biohit®, BRAND®, Capp®, Eppendorf®, Gilson®, Hamilton®, HTL, Nichiryo®, Rainin®, Sartorius®, Socorex®, STARLAB®, Thermo Scientific™