



# ErgoOne® E Single-Channel Pipette Electronic, 50 – 1000µl

SKU G9001-1000



## Product Highlights

- Reliable
- Lightweight, ergonomic design
- Adjustable finger rest for individual comfort
- Significantly reduced tip loading and ejection forces
- Universal tip fit
- Easy to use, intuitive menu structure
- Autoclavable pipette shaft
- 1-year warranty

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**STARLAB International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Channel type	1-channel
Color code	Blue
Volume Setting	Variable
Volume	1000 µl
Volume range	50 – 1,000 µl
Available options	12-channel, 8-channel, 1-channel, Multi-channel, Single-channel
Pack Size	1 Box × 1 Piece

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**STARLAB International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)




## More informations about ErgoOne® E Single-Channel Pipette Electronic, 50 – 1000µl









### Reliable and ergonomic!

ErgoOne® E single-channel pipettes provide reliable pipetting every day. Ergonomic and lightweight, ErgoOne® E pipettes are perfect for comfortable pipetting during long periods of use.

Five operating modes provide standard and reverse pipetting, mixing and dispensing, plus a unique “Electrophoresis” mode to help you load gels easily. The long lasting battery can also be charged during use so you are never without your pipette.



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Individual Stand (for 5000 µl Single-Channel Pipette)</a>	1 Box × 1 Piece	G9999-1003
 <a href="#">Charging Stand with AC Adapter</a>	1 Box × 1 Piece	G9999-1001
 <a href="#">Individual Stand (for Single-Channel Pipettes, 10 µl to 1000 µl)</a>	1 Box × 1 Piece	G9999-1002
 <a href="#">Fluorstatic Grease for Guiding Axis (for Multi-Channel Pipettes)</a>	1 Pack × 1 Piece	G9999-1008
 <a href="#">Silicone Oil for Piston (for 5000 µl Single-Channel and all Multi-Channel Pipettes)</a>	1 Pack × 1 Piece	G9999-1007
 <a href="#">Silicone Oil for Piston (for Single-Channel Pipettes, 10 µl to 1000 µl)</a>	1 Pack × 1 Piece	G9999-1006
 <a href="#">Replacement Battery for ErgoOne® E</a>	1 Box × 1 Piece	G9999-1005
 <a href="#">AC Adapter for Europe/UK (230V/50 Hz)</a>	1 Box × 1 Piece	G9999-1000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**STARLAB International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)