



## Plate Centrifuge

SKU N2631-0008



### Product Highlights

- Less than quarter the size of most plate centrifuges
- Accommodates two plates
- Rapid start to stop
- Use before and after thermal cycling to increase PCR yield
- Accepts skirted, semi-skirted and non-skirted plates
- Easy cleaning and maintenance
- 2-year warranty

## Product Applications

Quick and easy centrifugation of samples in PCR plates or microplates

Quick spin down of droplets and condensation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**STARLAB International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Pack Size	Box
Weight	1.5 kg
Rotor capacity	Two plates, any style, up to 25 mm high
Dimensions (W × L)	190 mm ×
Power supply	230 V, 50/60 Hz
Fuse	0.5 AT
Permissible Relative Humidity	< 85%
Dimensions (W × H)	190 × 193 mm
Permissible ambient temperature	4 to 30 °C
RCF	500 × g
Packaging Item	Piece
Packaging Unit	Box
Speed	2,500 rpm
Pack Size	1 Box × 1 Piece
Dimensions (W × D × H)	190 mm × 210 mm × 193 mm
External width	190 mm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**STARLAB International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)




## More informations about Plate Centrifuge

### You spin me round!

The plate centrifuge makes the safe and quick spin down of droplets and condensation very easy. Sealed plates are loaded vertically and surface tension keeps samples in place so there is no danger of spilling samples. The centrifuge accepts all plates up to 25 mm high. An optional adaptor is available for centrifuging 48-well PCR plates and 0.2 ml PCR tubes or strips.



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Adapter for 1/2 plates, strips and tube</a>	1 Box × 2 Piece	N2631-0018

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**STARLAB International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)